









Product Code . SEL-TM-11061

Drop-Point Apparatus

Description

Drop-Point Apparatus

The apparatus consists of brass sleeve and case with metal cup and a glass boiling tube with cork fitted to a bath (Beaker) is provided a stirrer is also provided in the bath. It is an indication of the critical temperature of the solid or semi solid to the liquid state.

We are leading manufacturers, suppliers of Drop-Point Apparatus for Thermodynamics Models. Contact us to get high quality Drop-Point Apparatus for Thermodynamics Models for schools, colleges, universities, research labs, laboratories and various industries.

{ "@context": "https://schema.org/", "@type": "Product", "name": "Drop-Point Apparatus", "image": "htt p://www.schooleducationalinstrument.com/images/catalog/product/2016120730Drop-PointApparatusWithlogo.jpg", "description": "The apparatus consists of brass sleeve and case with metal cup and a glass boiling tube with cork fitted to a bath (Beaker) is provided a stirrer is also provided in the bath. It is an indication of the critical temperature of the solid or semi solid to the liquid state. We are leading manufacturers, suppliers of Drop-Point Apparatus for Thermodynamics Models. Contact us to get high quality Drop-Point Apparatus for Thermodynamics Models for schools, colleges, universities, research labs, laboratories and various industries.", "brand": "School Lab Instrument", "sku": "5", "gtin8": "5", "gtin13": "5", "gtin14": "5", "mpn": "5", "aggregateRating": "0", "ratingCount": "15" } }

Direct Contact Details \(\subseteq +91-8569909696 \) sales@schooleducationalinstrument.com www.schooleducationalinstrument.com