



Product Code . SEL-LG-11278

Adapter Splash Head

Description

Adapter Splash Head

Prevents contents of the flask from being drawn into condenser.




For use with rotary evaporators.

Modified bump traps provide expanded volume for reactions that have tendency of excessive foam.

Material : Borosilicate Glass.

We are leading manufacturers, suppliers of Adapter Splash Head for Laboratory Glassware. Contact us to get high quality Adapter Splash Head for Laboratory Glassware for schools, colleges, universities, research labs, laboratories and various industries.

```
{ "@context": "https://schema.org/", "@type": "Product", "name": "Adapter Splash Head", "image": "http://www.schooleducationalinstrument.com/images/catalog/product/494701977AdapterSplashHead.jpg", "description": "Adapter Splash Head Prevents contents of the flask from being drawn into condenser. For use with rotary evaporators. Modified bump traps provide expanded volume for reactions that have tendency of excessive foam. Material : Borosilicate Glass. We are leading manufacturers, suppliers of Adapter Splash Head for Laboratory Glassware. Contact us to get high quality Adapter Splash Head for Laboratory Glassware for schools, colleges, universities, research labs, laboratories and various industries.", "brand": "School Lab Instrument", "sku": "5", "gtin8": "5", "gtin13": "5", "gtin14": "5", "mpn": "5", "aggregateRating": { "@type": "AggregateRating", "ratingValue": "5", "bestRating": "5", "worstRating": "0", "ratingCount": "15" } }
```

School Educational Instruments,
Hargolal Road, Ambala Cantt, Haryana India
Direct Contact Details  +91-8569909696  sales@schooleducationalinstrument.com
 www.schooleducationalinstrument.com