

Description

Chromatography Assembly, Consisting of Dropping Funnel

Reduction adapter with 19/26 socket
and 24/29 cone.

Column tube 2 cm. dia 20 cm.

Effective length with 24/29 socket, 19/26
cone and fused sintered disc.

Chromatography Assembly, "Lab
Glassware", consisting of dropping
funnel 250 ml. with 19/26 cone.

Adapter with tap 19/26 socket 24/29
cone.

Filtering flask 250 ml with 24/29 socket.

Tamping rod 50 cm. Long.

We are leading manufacturers, suppliers
of Chromatography Assembly,
Consisting of Dropping Funnel for
Laboratory Glassware. Contact us to get
high quality Chromatography Assembly,
Consisting of Dropping Funnel for
Laboratory Glassware for schools,
colleges, universities, research labs,
laboratories and various industries.

School Educational Instruments,
Hargolal Road, Ambala Cantt, Haryana
India

Direct Contact Details 📞

+91-8569909696 ✉ sales@schooleduca
tionalinstrument.com



www.schooleducationalinstrument.com